

Personalised Healthcare (PHC) with Foundation Medicine (FMI)

Fatma Elçin KINIKLI,

FMI Turkey, Science Leader

Agenda

- 
- **PHC Approach Provides Better Patient Outcome**

- 
- **FMI offers Comprehensive Genomic Profiling, providing patients with more and personalised therapeutic options to improve outcomes**

Agenda

- **PHC Approach Provides Better Patient Outcome**


- **FMI offers Comprehensive Genomic Profiling, providing patients with more and personalised therapeutic options to improve outcomes**

The genetic complexity of cancer

Cancer care becoming more personalised by understanding disease biology

Advances in cancer biology & diagnostics drive future options for cancer care


One drug fits all



One patient group,
one biomarker, one drug




One patient, one tumor
profile, more options




1. Agyeman, A.A. and Ofori-Asenso, R. (2015) *J Pharm Bioallied Sci* 7(3):239–244;
2. Schwaederle, M. and Kurzrock, R. (2015) *Oncoscience* 2(10):779-80.

Cancer care is increasingly complex



2000




2015

The classical approach to cancer is evolving




2000: Breast cancer subtypes



Clinical decisions based on **affected tissue, histology and disease stage**

2017: Genomic drivers




Clinical decisions based on the **results of comprehensive genomic profiling**

Unmet need in lung cancer

Negligible difference between chemotherapy regimens

Overall Survival (%)


Advanced NSCLC



Schiller JH et al. (2002) NEJM 346:92-98


Profiling has been shown to improve outcomes in lung cancer patients

Metastatic lung adenocarcinoma patients



- **Biomarkers associated with targeted therapies in guidelines:** HER2, EGFR, ALK, BRAF, MET
- **Biomarkers not associated with targeted therapies in guidelines:** PIK3CA, MEK, AKT1, KRAS, NRAS

Advanced NSCLC patients



- **Biomarkers associated with targeted therapies in guidelines:** HER2, EGFR, ALK, BRAF
- **Biomarkers not associated with targeted therapies in guidelines:** PIK3CA, KRAS

Meta analysis of ~85,000 patients in Phase I, II and III trials



Targeted therapy in of itself is not generally effective – a biomarker is needed

Worst outcome



Best outcome

ARMS type	Pooled Analysis			Meta-analysis		
	RR (%)	PFS (Mos)	OS (Mos)	RR (%)	PFS (Mos)	OS (Mos)
Non-personalized targeted	4	2.6	8.7	7.5	2.5	8.3
Cytotoxic	12	3.3	9.4	16.1	3.3	9.3
Personalized targeted	30	6.9	15.9	31.3	6.1	13.7

Agenda

- **PHC Approach Provides Better Patient Outcome**

- **FMI** offers Comprehensive Genomic Profiling, providing patients with more and personalised therapeutic options to improve outcomes

Roche and Foundation Medicine Inc. (FMI)

Entered into a broad strategic collaboration in 2015



The strategic collaboration aims to further advance **FMI's leading position in molecular information** and **analysis** while **providing Roche** a unique **opportunity to identify and develop novel treatment options** for **patients**

Foundation Medicine is a **molecular information company leading a transformation in cancer care**, where each patient's treatment is informed by a **deep understanding** of the **molecular changes** that contribute to their disease.

Foundation Medicine's comprehensive genomic profiling approach



**Data
aggregation
and analysis**



**Scientific/clinical
expert review**

FMi is taking patients from comprehensive identification of gene alterations to more personalised therapies and improved outcomes

The screenshot shows a patient report from Foundation Medicine. It includes patient information, a table of genomic alterations, and a table of therapeutic implications.

Gene	Alteration	Frequency	Therapeutic Implications
BRCA1	Pathogenic	100%	Yes, see clinical trials
BRCA2	Pathogenic	100%	Yes, see clinical trials
PTEN	Pathogenic	100%	Yes, see clinical trials
APC	Pathogenic	100%	Yes, see clinical trials
TP53	Pathogenic	100%	Yes, see clinical trials
MDM2	Pathogenic	100%	Yes, see clinical trials


**A report connecting
patients to targeted
therapies**

FMi: Foundation Medicine, Inc.; NGS: next-generation sequencing.

Foundation Medicine, Inc. (2017) <https://www.foundationmedicine.com/> Accessed Feb 2017; Foundation Medicine, Inc. Patient report.

First page of the FoundationOne® Report

MSI & TMB



Patient Name		Report Date	Tumor Type Skin melanoma
Date of Birth	Sex	Medical Facility	Specimen Received
	Male	Ordering Physician	Specimen Site
FMI Case #		Additional Recipient	Skin
Medical Record #		Medical Facility ID #	Date of Collection
Specimen ID		Pathologist	Specimen Type

ABOUT THE TEST:
FoundationOne™ is a next-generation sequencing (NGS) based assay that identifies genomic alterations within hundreds of cancer-related genes.

PATIENT RESULTS

- 10 genomic findings
- 10 therapies associated with potential clinical benefit
- 0 therapies associated with lack of response
- 20 clinical trials

TUMOR TYPE: SKIN MELANOMA

Genomic Alterations Identified[†]

NF1 splice site 3709-2A>G
NRAS Q61K
MET R1148Q
PTEN H93Y
CDKN2A p14ARF I27fs*20
KDM5C S539*
KEL A17E – subclonal*
RB1 K202*
TERT promoter -124C>T

Additional Findings[†]

Tumor Mutation Burden TMB-High; 34.34 Muts/Mb

Additional Disease-relevant Genes with No Reportable Alterations Identified[†]

BRAF
KIT

[†] For a complete list of the genes assayed and performance specifications, please refer to the Appendix
^{*} See Appendix for details

THERAPEUTIC IMPLICATIONS

Genomic Findings Detected	FDA-Approved Therapies (in patient's tumor type)	FDA-Approved Therapies (in another tumor type)	Potential Clinical Trials
<i>Tumor Mutation Burden</i> TMB-High; 34.34 Muts/Mb	Pembrolizumab Nivolumab Ipilimumab	Atezolizumab	Yes, see clinical trials section
<i>NF1</i> splice site 3709-2A>G	Cobimetinib Trametinib	None	Yes, see clinical trials section
<i>NRAS</i> Q61K	Cobimetinib Trametinib	None	Yes, see clinical trials section

For more comprehensive information please log on to the Interactive Cancer Explorer™
 To set up your Interactive Cancer Explorer account, contact your sales representative or call 1-888-988-3639.



FOUNDATION
MEDICINE®



Section on identified genomic alterations

The **GENOMIC ALTERATIONS** section contains a detailed interpretive statement that references:

- Each **gene** and its alteration(s)
- **Frequency** of each genomic alteration(s) identified in patient cohorts with similar tumor type
- **Prognosis** and biological implication of each genomic alteration(s)
- Potential **treatment strategies** specific for each genomic alteration

**FOUNDATION**
ONE

Patient Name
Not Given

Report Date
30 August 2014

Tumor Type
Lung
adenocarcinoma

GENOMIC ALTERATIONS

GENE ALTERATION	INTERPRETATION
MET • exon 14 splice site (3028+2T>A)	<p>Gene and Alteration: MET, also known as hepatocyte growth factor receptor (HGFR) or c-MET, encodes a receptor tyrosine kinase that is activated by the ligand HGF. MET activation results in signaling mediated in part through the RAS-RAF-MEK and PI3K pathways to promote proliferation¹. MET alterations predicted to affect exon 14 splicing, such as seen here, have been reported in 1.4% of lung tumors^{2,3,4,5,6} and have been associated with skipping of exon 14^{7,8,9}. Loss of exon 14 increases MET stability, leading to prolonged signaling upon HGF stimulation and increased oncogenic potential^{10,11,12,13}. This mutation is therefore expected to be activating. Responses to various MET inhibitors have been reported for multiple patients with alterations in their tumors predicted to affect splicing of MET exon 14^{2,3,4,14,15}.</p> <p>Frequency and Prognosis: In one study of 4432 lung adenocarcinoma cases, MET mutations (primarily those affecting MET exon 14 splicing) have been reported in ~3% of samples². In the TCGA datasets, MET mutation has been observed in 8.3% of lung adenocarcinomas and 2.1% of lung squamous cell carcinomas¹⁷. MET amplification has been reported at incidences of 14-48% in NSCLC and correlated with increased MET protein expression^{18,19,20,21,22}. Studies into the effect of MET amplification on prognosis for patients with NSCLC have yielded conflicting results^{18,22,23,24,25}. One study observed an association between MET amplification and increased long-term survival for patients with lung adenocarcinoma, although concurrent MET amplification and EGFR mutation have been correlated with reduced disease-free survival^{27,28}.</p> <p>Potential Treatment Strategies: MET amplification or activating mutations may predict sensitivity to targeted therapies²⁹, such as the kinase inhibitors crizotinib and cabozantinib. Crizotinib is FDA approved for the treatment of ALK-positive NSCLC³⁰. Cabozantinib is FDA approved for the treatment of metastatic medullary thyroid cancer³¹, and a patient with a MET amplification and a mutation associated with MET exon 14 splicing has achieved a complete response to cabozantinib³². Crizotinib has benefited patients with MET-amplified NSCLC³³, lung squamous cell carcinoma³⁴, lung adenocarcinoma³⁵, gastroesophageal cancer³⁶, glioblastoma³⁷, and carcinoma of unknown primary³⁷. In one clinical trial, treatment with AMG 337 led to a response rate of 50% (5/10), including 1 complete response, for patients with MET-amplified gastric, esophageal, or gastroesophageal junction cancer (Kwak et al., 2015, ASCO GI Abstract 01). MET-targeting antibodies onartuzumab and MetMab have elicited responses in patients with MET-amplified NSCLC³⁸ and gastric cancer³⁹. In addition, high MET expression has been suggested to predict patient response to therapy regimens involving rituximab, a monoclonal HGF-targeting antibody (Oliner et al., 2012, ASCO Abstract 4005). Furthermore, MET inhibitors crizotinib, capmatinib, PF-04217903, and foretinib have provided benefit to patients with MET-mutated papillary renal cell carcinoma^{40,41}, histiocytic sarcoma⁴², lung adenocarcinoma^{3,14,15}, lung large cell carcinoma⁴³, and lung squamous cell carcinoma⁴⁴.</p>
CDKN2A/B • loss	<p>Gene and Alteration: CDKN2A encodes two different, unrelated tumor suppressor proteins, p16INK4a and p14ARF, whereas CDKN2B encodes the tumor suppressor p15INK4b^{45,46}. Both p15INK4b and p16INK4a bind to and inhibit CDK4 and CDK6, maintaining the growth-suppressive activity of the RB tumor suppressor; loss or inactivation of either p15INK4b or p16INK4a contributes to dysregulation of the CDK4/6-cyclin D-RB pathway and loss of cell cycle control^{46,47}. The p14ARF tumor-suppressive functions</p>

For more comprehensive information please log on to the Interactive Cancer ExplorerTM
To set up your Interactive Cancer Explorer account, contact your sales representative or call (888) 988-3839.

Electronically Signed by | CJA Number: 2020270702 | 28 August 2014
Foundation Medicine, Inc., 150 2nd Street, 1st Floor, Cambridge, MA 02141 | 1.800.898.3839

page 3 of 21



FOUNDATION
MEDICINE



THERAPEUTIC IMPLICATIONS

Genomic Findings Detected	FDA-Approved Therapies (in patient's tumor type)	FDA-Approved Therapies (in another tumor type)	Potential Clinical Trials
<i>CCND1</i> amplification	None	Palbociclib Ribociclib	Yes, see clinical trials section
<i>PIK3CB</i> amplification	None	None	Yes, see clinical trials section
<i>EMSY</i> amplification	None	None	None
<i>FGF19</i> amplification	None	None	
<i>FGF3</i> amplification	None	None	

 **PIK3CB**
amplification

Gene and Alteration: PIK3CB (PI3K-beta, p110-beta) encodes a member of the class IA phosphoinositide 3-kinases (PIK3CA, PIK3CB, PIK3CD, and PIK3CG), which are essential regulators of cellular proliferation, survival, metabolism and motility. Dysregulation of the PI3K signaling pathway is observed in many human cancers and occurs most frequently through loss of the tumor suppressor PTEN or activating mutations in PIK3CA (p110-alpha). Contrary to PIK3CA, which is frequently mutated within hot spots causing increased kinase activity, mutation of PIK3CB is quite rare^{20,21,22}. PIK3CB amplification has been identified in various tumor types^{9,23,24}, and PIK3CB overexpression has been shown to be oncogenic in vitro²⁵. A patient with PTEN-deficient prostate cancer and PIK3CB amplification had a partial response to the PI3K-beta-selective inhibitor GSK2636771 (de Bono et al., 2015; AACR Abstract CT328).

Frequency and Prognosis: PIK3CB mutation has been reported in up to 1% of prostate adenocarcinoma cases, while amplification has been found in up to 2% of samples analyzed^{5,26}. A constitutively active form of p110-beta induces prostatic intraepithelial neoplasia in a murine prostate tumor model²⁷. In addition, inactivation of p110-beta, but not p110-alpha, was shown to inhibit prostate cancer development in a mouse model drive by PTEN heterozygous deletion²⁸.

Potential Treatment Strategies: Preclinical evidence indicates that PTEN-deficient tumors depend on PI3K-beta^{29,30,31}, and inhibitors selective for PI3K-beta, including GSK2636771 and AZD8186, are in clinical trials for PTEN-deficient tumors (Arkenau et al., 2014; ASCO Abstract 2514, Siu et al., 2015; AACR Abstract CT329). PIK3CB activating alterations may further refine the use of PI3K-beta-selective inhibitors²². Whereas two patients with PTEN-deficient prostate tumors and concurrent PIK3CB L1049R mutation or PIK3CB amplification derived significant clinical benefit with GSK2636771 (de Bono et al., 2015; AACR Abstract CT328), PIK3CB D1067 mutations rendered PTEN-deficient breast cancer cells less sensitive to pan-PI3K or PI3K-beta-selective inhibitors in preclinical assays³².

Section on possible treatments

The THERAPIES section provides further details on:

- **Approved therapies** to which the patient's tumor type may be sensitive or resistant based on genomic profile
- Approved therapies associated with benefit based on similar genomic alterations in **other tumor types***

With information on approved indications, gene association and supporting data

* IMPORTANT: this does **not** indicate evidence for off-label use in patient's tumor type

FOUNDATION ONE		Patient Name Not Given	Report Date 20 August 2014	Tumor Type Lung adenocarcinoma
THERAPIES				
FDA APPROVED THERAPIES IN PATIENT TUMOR TYPE				
THERAPY SUMMARY OF DATA IN PATIENT TUMOR TYPE				
Crizotinib				
<p>Approved Indications: Crizotinib is an inhibitor of the kinases MET, ALK, ROS1, and RON. It is FDA approved to treat non-small cell lung cancer (NSCLC) in patients whose tumors are positive for ALK rearrangements as detected by an FDA-approved test.</p> <p>Gene Association: Activating MET alterations may confer sensitivity to crizotinib. Crizotinib has benefited patients with NSCLC or histiocytic sarcoma tumors harboring various alterations associated with MET exon 14 skipping^{2,3,14,15,16}.</p> <p>Supporting Data: Crizotinib has demonstrated efficacy in patients with NSCLC and ALK rearrangements¹⁰, ROS1 rearrangements¹⁹, or MET activation^{2,3,14,15,16,20}. Crizotinib has benefited patients with MET-amplified NSCLC²¹, lung squamous cell carcinoma²², and lung adenocarcinoma²³ and achieved partial responses for three patients with lung adenocarcinoma and MET exon 14 splice site alterations²⁴.</p>				
ADDITIONAL THERAPIES – FDA APPROVED IN OTHER TUMOR TYPES				
THERAPY SUMMARY OF DATA IN OTHER TUMOR TYPE				
Cabozantinib				
<p>Approved Indications: Cabozantinib is a kinase inhibitor that targets MET, RET, VEGFRs, KIT, FLT-3, TIE-2, AXL, and TRKB. It is FDA approved for the treatment of medullary thyroid cancer.</p> <p>Gene Association: Alterations that are expected to increase signaling through these kinases, in particular MET, RET, and VEGFR-2, may predict sensitivity to treatment with cabozantinib. A patient with a MET amplification and a mutation associated with MET exon 14 skipping has achieved a complete response to cabozantinib⁷.</p> <p>Supporting Data: A preclinical study in mouse models of multiple tumor types, including lung cancer, reported that cabozantinib treatment resulted in decreased cell proliferation and tumor growth as well as increased apoptosis¹⁰. A Phase 2 randomized discontinuation trial of cabozantinib in metastatic non-small cell lung cancer (NSCLC) reported a 54% rate of tumor regression in heavily pretreated patients, with a safety profile similar to that of other tyrosine kinase inhibitors (Hellestedt et al., 2012; ASCO Abstract 7514). Preclinical studies have shown that lung cancer cells with dual EGFR and MET activation respond to combination treatment with EGFR inhibitors and MET inhibitors^{18,19}, and clinical trials investigating this combination are under way in NSCLC. A Phase 2 randomized trial of cabozantinib, erlotinib, or a combination in patients with EGFR-wild-type NSCLC reported improved PFS (3.9 months vs. 1.9 months on erlotinib alone, HR 0.33, p = 0.0002) and OS (HR 0.50, p = 0.02) for cabozantinib treatment compared to erlotinib, although MET expression did not correlate with response (Neel et al., 2015; ASCO Abstract 8003). As part of an ongoing Phase 2 study of cabozantinib for the treatment of NSCLC, 3 partial responses and two instances of stable disease have been achieved by patients with RET rearrangements in their tumors^{25,26}.</p>				
For more comprehensive information please log on to the Interactive Cancer Explorer™ To set up your Interactive Cancer Explorer account, contact your sales representative or call (888) 988-3639. Devisively Signed by CLIA Number: 2202827021 20 August 2014 Foundation Medicine, Inc., 150 2nd Street, 1st Floor, Cambridge, MA 02141 1.888.988.3639 page 6 of 21				

Section on clinical trials to consider

The CLINICAL TRIALS TO CONSIDER section of the report provides:

- **Rationale** for potential clinical trials
- Details on **keyword terms** used for search on [clinicaltrials.gov](https://www.clinicaltrials.gov) for relevant clinical trials
- Details on currently available clinical trials for which the patient **may be eligible** based on the genomic profile of patient's tumor

Patient Name
Not Given

Report Date
20 August 2014

Tumor Type
Lung
adenocarcinoma

CLINICAL TRIALS TO CONSIDER

IMPORTANT: While every effort is made to ensure the accuracy of the information contained below, the information available in the public domain is continuously updated and should be investigated by the physician or research staff. This is not meant to be a complete list of available trials. In order to conduct a more thorough search, please go to www.clinicaltrials.gov and use the search terms provided below. For more information about a specific clinical trial, type the NCT ID of the trial indicated below into the search bar.

GENE

RATIONALE FOR POTENTIAL CLINICAL TRIALS

MET
exon 14 splice site (2028:27:TA)

MET amplification or activating mutations may confer sensitivity to therapies targeting this kinase. Examples of clinical trials that may be appropriate for this patient are listed below. These trials were identified through a search of the trial website [clinicaltrials.gov](https://www.clinicaltrials.gov) using keyword terms such as "MET", "crizotinib", "vandetanib", "NSCLC", "lung", "solid tumor", and/or "advanced cancer".

TITLE	PHASE	TARGETS	LOCATIONS	NCT ID
Phase 1 Safety, Pharmacokinetic And Pharmacodynamic Study Of PF-02341066, A o-Met/KIT/FGFR Selective Tyrosine Kinase Inhibitor, Administered Orally To Patients With Advanced Cancer	Phase 1	MET	California, Colorado, Illinois, Massachusetts, Michigan, New York, North Carolina, Ohio, Pennsylvania, Tennessee, Seoul (Korea, Republic of), Victoria (Australia)	NCT002585195
Phase I Study of INC280 Plus Erlotinib in Patients With C-Met Expressing Non-Small Cell Lung Cancer	Phase 1	MET, EGFR	California	NCT01811307
A Phase I Open-label Dose Escalation Study With Expansion to Assess the Safety and Tolerability of INC280 in Patients With o-MET Dependent Advanced Solid Tumors	Phase 1	MET	Arkansas, Illinois, Maryland, Michigan, New York, Tennessee, Texas, Utah, multiple ex-US locations	NCT01324470
Open-Label Dose-Escalation Trial to Evaluate the Safety, Pharmacokinetics, and Pharmacodynamics of Daily Oral MGC0295 Administered Without Interruption to Subjects With Advanced Malignancies	Phase 1	AXL, DDR2, Eph, FLT1, IGF1R, KIT, MET, PDGFRb, RET, TRK, VEGFRs	California, Illinois, Massachusetts, Missouri, New York, North Carolina, Pennsylvania, Texas, Utah, Washington, Alberta (Canada), British Columbia (Canada), Gyeonggi-do (Korea, Republic of), Quebec (Canada), Seoul (Korea, Republic of)	NCT00967832
A Phase II, Multicenter, Open-label Study of EGFR16 in Combination With Nivolumab in Adult Patients With EGFR Mutated Non-small Cell Lung Cancer and of INC280 in Combination With Nivolumab in Adult Patients With oMet Positive Non-small Cell Lung Cancer	Phase 2	MET, EGFR, PD-1	Amsterdam (Netherlands), Chur (Switzerland), New South Wales (Australia), PG (Italy), PI (Italy), PN (Italy), Queensland (Australia), Singapore (Singapore), South Australia (Australia)	NCT02323126

For more comprehensive information please log on to the Interactive Cancer Explorer™
To set up your Interactive Cancer Explorer account, contact your sales representative or call (888) 988-3639.

Devisively Signed by | CLIA Number: 2202027021 | 20 August 2014
Foundation Medicine, Inc., 150 2nd Street, 14th Floor, Cambridge, MA 02141 | 1.888.988.3639

page 8 of 21



**FOUNDATION
MEDICINE®**




Agenda

- **PHC Approach Provides Better Patient Outcome**

- **FMI offers Comprehensive Genomic Profiling,** providing patients with more and personalised therapeutic options to improve outcomes

Four types of ways genes can be altered




Routine single-marker molecular test

The most common type of molecular testing

CATEGORY ONE

Routine single marker molecular tests such as **IHC**, **PCR** and **FISH** that have been used for decades and will continue to play an important

CATEGORY TWO
on cancer diagnosis
CATEGORY THREE



FISH: fluorescence *in situ* hybridization; IHC: Immunohistochemistry; PCR: Polymerase chain reaction



FOUNDATION
MEDICINE®



Multi-gene “hot spot” test

Broader testing focused on a narrow subset of

genes

CATEGORY ONE

CATEGORY TWO

The **hot spot**

NGS panels

identify pre-

specified

mutations

occurring in

very limited

areas of genes

of interest and

fail to detect all

CATEGORY THREE



NGS: Next-generation sequencing

Foundation Medicine®

The most comprehensive genomic test available

CATEGORY ONE

CATEGORY


TWO

CATEGORY THREE

FMI's

**comprehensive
genomic
profiling**

approach of
testing all of the
known clinically
relevant cancer
genes for **all**



FMI: Foundation Medicine, Inc.



FOUNDATION
MEDICINE®




Two Different Targeted Sequencing Approaches

Comprehensive Genomic Profiling detects more than Hot Spot NGS

Example:
EGFR
gene

Hot spot tests detect selective alterations in selective parts of the EGFR gene*



Agenda

1. Foundation Medicine, Inc. (2017) <https://www.foundationmedicine.com/> Accessed Feb 2017; 2. Schwaederle, M., et al. (2015) *Mol Cancer Ther.* 14(6):1488-94; 3. Wheler, J., et al. (2016) *Cancer Research* 76(13): 3690-701; 4. Rozenblum, A.B., et al. (2017) *J Thorac Oncol.* 12(2):258-68; 5. Drilon, A., et al. (2015) *Clin Cancer Res.* 21(16):3631-9. 6. Kris MG et al. (2014) *JAMA* 311(10):1098-2006. 7. Barlesi F et al. (2016) *Lancet* S0140-6736(16)

MSKCC experience with FoundationOne®


Supporting first-line use of CGP in lung cancer

Background

- **31** lung adenocarcinoma patients with available tissue specimen were identified
- Those patient's tumors harbored **no evidence of a genomic alteration** via extensive, focused non-NGS testing covering known lung cancer alterations in 11 genes (*EGFR*, *ERBB2*, *KRAS*, *NRAS*, *BRAF*, *MAP2K1*, *PIK3CA*, *AKT1*, *ALK*)

Results

- FoundationOne® detected genomic alterations in **94% of patients with pan-negative lung adenocarcinoma**



These findings support the **first-line use** of CGP to detect the broadest range of genomic alterations in lung adenocarcinomas while **preserving tissue**

GA: Genomic Alteration

Clinical utility of finding more alterations with FoundationOne[®]

NSCLC patients can benefit from targeted therapies

1
7
%

of EGFR exon 19 deletions
missed by hotspot tests¹

3
5
%

of ALK-rearranged cases
missed by FISH³

7
5
%

of NSCLC patients with EGFR
exon 19 deletions can
respond to EGFR tyrosine
kinase inhibitors, with median
OS > 1 year²

8
0
%

of ALK-rearranged patients
identified by FoundationOne
respond to ALK inhibitor
crizotinib³

1. Schrock AB et al. (2016) *Clin Cancer Res*. Mar 1.
2. Sequist LV et al. (2007) *J Clin Oncol*. 25:587-95.
3. Ali AM et al. (2016) *The Oncologist*

Patient Case

Patient Information



- 58-year-old female

Diagnosis

- Triple negative inflammatory breast cancer


Treatment status

- Refractory to taxanes, dasatinib, bevacizumab, ixabepilone and gemcitabine
- After progressive disease additional biomarker testing done and histology, ER/PR (by IHC) and HER2 (by FISH) confirmed as all negative



FoundationOne® analysis

- The same sample sent for CGP. HER2 was consistent with FISH, but **two distinct HER2 mutations** (V777L and S310F) found



Patient Case

Patient Information



- 43-year-old female

Diagnosis

- **Breast carcinoma.**
- **ER** positive, **PR** negative and **HER2 negative** by **FISH** and **IHC**.

Treatment status

After chemotherapy peritoneal and omental metastases, a left adrenal mass and enlargement of liver metastases are observed.

FoundationOne® analysis

The sample was found to have an ERBB2 L755S base substitution mutation and the tumors responded to neratinib

TABLE 2. Summary of Clinical Characteristics for 138 Breast Cancers With ERBB2 Short-Variant Mutations

	NOS	IDC	ILC	MUC	
No. of cases	69	40	27	2	
Median age	61 y	60 y	58 y	58 y	
Average age	60.9 y	58.5 y	58 y	58 y	
Specimen site					Percentage of all ERBB2 mut samples
Bone	0	1	0	0	0.7%
Breast	7	28	16	0	37%
Liver	18	5	3	0	18.8%
Lung	5	0	1	0	4.3%
Lymph node	14	2	3	1	14.5%
Other	24	5	4	1	24.6%

Patient Case

Patient Information



- 58-year-old female

Diagnosis

- **Metastatic inflammatory breast carcinoma.**
- **Triple negative** by **FISH** and **IHC**.

Treatment status

Patient was previously treated with multiple courses of cytotoxic chemotherapy. However disease remain active.

FoundationOne® analysis

2 distinct ERBB2mut, 1 each in the kinase domain (V777L) and the ECD (S310F). The patient received multiple anti-HER2 targeted therapies combined with chemotherapy and demonstrated substantial clinical response.

Patient Case

Patient Information



- 58-year-old female

Diagnosis

- **Breast carcinoma.**
- **ER** positive, **PR** positive and **HER2 negative** by **FISH** and **IHC**.
- With the treatment refractory hepatic metastasis.

FoundationOne® analysis

The sample was found to have ERBB2 S310F mutation and positive HER2.




Figure 2. Long tail plot of genomic alterations in 138 cases of ERBB2-mutated recurrent metastatic breast cancer.

Agenda

1. Foundation Medicine, Inc. (2017) <https://www.foundationmedicine.com/> Accessed Feb 2017; 2. Schwaederle, M., et al. (2015) *Mol Cancer Ther.* 14(6):1488-94; 3. Wheler, J., et al. (2016) *Cancer Research* 76(13): 3690-701; 4. Rozenblum, A.B., et al. (2017) *J Thorac Oncol.* 12(2):258-68; 5. Drilon, A., et al. (2015) *Clin Cancer Res.* 21(16):3631-9. 6. Kris MG et al. (2014) *JAMA* 311(19):1998-2006. 7. Barlesi F et al. (2016) *Lancet* S0140-6736(16)

Analytic Validation

nature
biotechnology

Development and validation of a clinical cancer genomic profiling test based on massively parallel DNA sequencing

Garrett M Frampton^{1,9}, Alex Fichtenholtz^{1,9}, Geoff A Otto¹, Kai Wang¹, Sean R Downing¹, Jie He¹, Michael Schnall-Levin¹, Jared White¹, Eric M Sanford¹, Peter An¹, James Sun¹, Frank Juhn¹, Kristina Brennan¹, Kiel Iwanik¹, Ashley Mailet¹, Jamie Buell¹, Emily White¹, Mandy Zhao¹, Sohail Balasubramanian¹, Selmira Terzic¹, Tina Richards¹, Vera Banning¹, Lazaro Garcia¹, Kristen Mahoney¹, Zac Zwirko¹, Amy Donahue¹, Himisha Beltran^{2,3}, Juan Miguel Mosquera^{3,4}, Mark A Rubin^{3,4}, Snjezana Dogan⁵, Cyrus V Hedvat⁵, Michael F Berger⁵, Lajos Pusztai⁶, Matthias Lechner⁷, Chris Boshoff⁷, Mirna Jarosz¹, Christine Vietz¹, Alex Parker¹, Vincent A Miller¹, Jeffrey S Ross^{1,8}, John Curran¹, Maureen T Cronin¹, Philip J Stephens¹, Doron Lipson¹ & Roman Yelensky¹

As more clinically relevant cancer genes are identified, comprehensive diagnostic approaches are needed to match patients to therapies, raising the challenge of optimization and analytical validation of assays that interrogate millions of bases of cancer genomes altered by multiple mechanisms. Here we describe a test based on massively parallel DNA sequencing to characterize base substitutions, short insertions and deletions (indels), copy number alterations and selected fusions across 287 cancer-related genes from routine formalin-fixed and paraffin-embedded (FFPE) clinical specimens. We implemented a practical validation strategy with reference samples of pooled cell lines that model key determinants of accuracy, including mutant allele frequency, indel length and amplitude of copy change. Test sensitivity achieved was 95–99% across alteration types, with high specificity (positive predictive value >99%). We confirmed accuracy using 249 FFPE cancer specimens characterized by established assays. Application of the test to 2,221 clinical cases revealed clinically actionable alterations in 76% of tumors, three times the number of actionable alterations detected by current diagnostic tests.

Base Substitutions

Sensitivity: >99.9% PPV: >99%

Insertions/Deletions

Sensitivity: 98% PPV: >99%

Copy Number Alterations

Sensitivity: >95% PPV: >99%

Gene Fusions¹

Sensitivity: >95% PPV: >99%
(>99% for ALK fusion²)

1. Based on analysis of coverage and re-arrangement structure in the COSMIC database for solid tumor fusion genes where alteration prevalence could be established, complemented by detection of exemplar rearrangements in cell line titration experiments.
2. Yelensky et al, *Presented at AACR 2014*



FOUNDATION
MEDICINE®



FoundationOne CDx : First Commercial Pan-tumour Comprehensive Genomic Profiling Assay Approved by FDA

- **First assay**, incorporating a broad range of **companion diagnostics**
- **17 targeted therapies** across five types of advanced cancers: non-small cell lung, melanoma, breast, colorectal and ovarian cancers
- Reports the genomic markers **MSI** and **TMB** to help inform decisions on **immunotherapy**



Summary

Why consider profiling with Foundation Medicine?

- Profiling has been shown to **improve outcomes for patients**
- **Foundation Medicine's** profiling services are designed to **capture all four types of genomic alterations** and accurately identifies actionable targets across a spectrum of cancers
- These alterations are delivered in a **comprehensive report** which describes potential therapies, trials, and the latest clinical literature to inform physician's decisions
- Evidence has shown **FoundationOne® detects alterations in patients that are pan-negative with single gene panels**, and in some indications can improve outcomes



Profiling with FoundationOne® finds more clinically-relevant alterations and can lead to better patient outcomes

1. Kris MG et al. (2014). JAMA 311(19):1998-2006; 2. Barlesi F et al. (2016). Lancet S0140-6736(16); 3. Ganesan P et al. (2014) Mol Cancer Ther; 13(12); 3175-84; 4. Ali S et al. (2016) Oncologist. 5 Schrock AB et al. (2016) Clin Cancer Res. Mar 1. pii: clincanres.1668.2015; 6 Suh J et al. (2016) Oncologist; 7. Drilon A et al. (2015) Clin Cancer Res 21(16):3631-9; 8.

***Doing now what patients need
next***